

Provide
sustainable
patency.



Resonance[®]

METALLIC URETERAL STENT

A clear introducer system

Resonance's stent positioner and clear sheath provide improved visualization and reference points for first and second pigtail deployment for all stent sizes.

A maximum indwelling time of 12 months

Resonance reduces the need for frequent stent changes and may thereby decrease the risk of infection.

MRI Conditional (up to 3 tesla)

Resonance has been demonstrated to pose no known hazard in a specified MRI environment with specified conditions of use.




www.cookmedical.com



Resonance[®]

METALLIC URETERAL STENT

The Resonance stent is specifically designed for extrinsic compression to mitigate many of the problems encountered with traditional plastic ureteral stents and other forms of treatment. The optimized compressive and radial strength of Resonance allows the stent to remain indwelling for up to 12 months. This long indwelling time reduces the need for frequent stent changes and may thereby decrease the risk of infection. The tightly wound metallic coil maintains its patency so that urine can drain continuously under severe extrinsic compression. Resonance is introduced coaxially through a sheath and removed using standard cystoscopic techniques.



Order Number	Reference Part Number	Fr	Length cm
Set			
G34125	RMS-060020	6.0	20
G34126	RMS-060022	6.0	22
G34127	RMS-060024	6.0	24
G34009	RMS-060026	6.0	26
G34128	RMS-060028	6.0	28
G34177	RMS-060030	6.0	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Order Number	Reference Part Number	Fr	Length cm
Stent and Positioning System Only			
G34108	RMS-060020-R	6.0	20
G34109	RMS-060022-R	6.0	22
G34110	RMS-060024-R	6.0	24
G34111	RMS-060026-R	6.0	26
G34112	RMS-060028-R	6.0	28
G34176	RMS-060030-R	6.0	30

Positioning System Only

G34178	STPV-083070-RMS
--------	-----------------

Customer Service - Urology

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

